

# MTXQCvX2 documentation

*Christin Zasada*

*2018-11-01*



# Contents

<b>1</b>	<b>Welcome</b>	<b>5</b>
<b>2</b>	<b>Introduction</b>	<b>7</b>
<b>3</b>	<b>Workflow for Maui-annotation projects</b>	<b>9</b>
3.1	Read this in case . . . . .	9
3.2	Quick view . . . . .	9
3.3	Input files . . . . .	9
3.4	Annotation-file . . . . .	10
3.5	Sample_extracts-file . . . . .	10
3.6	Internal Standard . . . . .	11
<b>4</b>	<b>Workflow for Metmax-extracted projects</b>	<b>13</b>
4.1	You want to follow this ... . . . .	13
4.2	Introduction . . . . .	13
4.3	Quick view . . . . .	13
4.4	Input files . . . . .	14
4.5	Annotation-file . . . . .	14
4.6	Sample_extracts-file . . . . .	14
4.7	Update metabolite names in <code>conversion_metabolite.csv</code> . . . . .	15
<b>5</b>	<b>MTXQCvX2 Universe</b>	<b>17</b>
5.1	MTXQC_init.Rmd . . . . .	17
5.2	MTXQC_ExperimentalSetup . . . . .	17
5.3	MTXQC_part 1 . . . . .	17
5.4	MTXQC_part 2 - Post-Processing . . . . .	17
5.5	MTXQC_part 3 - Manual Validation . . . . .	17
5.6	MTXQC_part 4 - Metmax integration . . . . .	17
<b>6</b>	<b>Configuration of MTXQCvX2 - <code>config_mtx/...</code></b>	<b>19</b>
6.1	<code>conversion_metabolite.csv</code> . . . . .	19
6.2	<code>letter_pathway_complete.csv</code> . . . . .	20
6.3	<code>quant1-values.csv</code> . . . . .	21
6.4	<code>incorporation_calc.csv</code> & <code>mid_backups.csv</code> . . . . .	21
<b>7</b>	<b>pSIRM experiments</b>	<b>23</b>
7.1	Experimental design . . . . .	23
7.2	Experimental procedures . . . . .	24
7.3	Protocol pSIRM . . . . .	25
7.4	Hints & notes . . . . .	26
<b>8</b>	<b>Experimental Protocols</b>	<b>27</b>
8.1	Metabolomics sample extraction . . . . .	27

8.2	Sample derivatisation . . . . .	29
<b>9</b>	<b>Data processing - MAUI</b>	<b>31</b>
9.1	Processing In ChromaToF . . . . .	31
9.2	Maui notes . . . . .	32
9.3	Maui exports . . . . .	32
<b>10</b>	<b>Data Processing - Metmax</b>	<b>33</b>
10.1	Resampling . . . . .	33
10.2	1D-basic . . . . .	33
10.3	Reference search . . . . .	33
10.4	Export for Metmax . . . . .	33
10.5	Data extraction with Metmax . . . . .	33
<b>11</b>	<b>Solutions and Materials</b>	<b>35</b>
11.1	Solutions . . . . .	35
11.2	Materials . . . . .	36
<b>12</b>	<b>Output tables</b>	<b>39</b>
12.1	Project-files . . . . .	39
12.2	output/gc/... . . . .	41
12.3	output/quant/...{o_quant} . . . . .	45
12.4	output/inc/... . . . .	56
<b>13</b>	<b>Frequently Asked Questions</b>	<b>59</b>
13.1	What are additional quantification standards . . . . .	59
13.2	How do I extend conversion_metabolite.csv . . . . .	59
<b>14</b>	<b>GC-MS settings</b>	<b>61</b>

# Chapter 1

## Welcome

This documentation introduced to you how to use MTXQCvX2 in order to assess the quality of your GC-MS derived data, perform the determination of calibration curves and absolute quantification. It furthermore provides you two normalisation strategies and the calculation of quantities in, e.g., pmol/1e+6 cells or pmol/mg tissue.

MTXQCvX2 does also enable the calculation of stable isotopic incorporation and the evaluation of the underlying data, the mass isotopomer distributions (MIDs).

The tool has been set up to support the in-lab developed workflow for quantitative metabolomics experiments using the in-house developed software Maui for the annotation of data. MTXQCvX2 bridges the gap between quality control and first data post-processing / analysis of GC-MS derived data (MTXQCvX2\_part1, MTXQCvX2\_part2).

Nevertheless MTXQCvX2 includes a module in order to integrate all kind of data provided in spreadsheet-format, e.g., derived from metmax, extracting required information and creating corresponding files (MTXQCvX2\_part4).

Both workflows are introduced in the distinct chapters including their required input parameter (chapter Workflow Maui and Workflow Metmax).



## Chapter 2

# Introduction

Experimental and mathematical concepts have been introduced for the pulsed stable isotope resolved metabolomics (pSIRM) approach in (Pietzke et al., 2014).

```
knitr::include_graphics("images/now.png")
```




Figure 2.1: Test external figure and how to include it



## Chapter 3

# Workflow for Maui-annotation projects

### 3.1 Read this in case

- you have run a Maui project
- exported all required container (see ??)
- you have a copy of sequence list and experimental conditions
- you know the extraction procedure

The following article describes briefly how to use MTXQCvX2 in case you used Maui for the annotation of your metabolomics project. It does not matter if you have performed an experiment including stable isotopes or if you just aim for the quantification of a few intermediates.

### 3.2 Quick view

1. Setup a new R-project and copy MTXQC template files and folders
2. Knit with parameter: `MTXQC_init.Rmd` and create project folder, e.g., `psirm_glucose`
3. Copy input files and rename `ManualQuantTable.tsv` (`e18205cz.tsv`)
4. Create `annotation.csv` and `sample_extracts.csv` files
5. Define the internal extraction standard
6. Knit with parameter: `MTXQC_ExperimentalSetup.Rmd`
7. Knit with parameter: `MTXQC_part1.Rmd`
8. Knit with parameter: `MTXQC_part2.Rmd`
9. If required, proceed with `MTXQC_part3.Rmd` for ManualValidation

### 3.3 Input files

Three different kind of export functions have been implemented in Maui. These functions provide the export of the actual data into `.csv` or `.tsv` files that are directly usable as input files for MTXQCvX2. Please refer to section 9.3 how you perform the export and which containers have to be exported using what export function and where to copy them in `psirm_glucose/input/`.

Certain circumstances might wish you to combine *multiple MAUI-projects* into one MTXQC-project. This might be the case when you run the same samples in split and splitless mode on the machine or your

experimental setup has been measured in multiple batches in order to avoid derivatisation effects.

It is highly recommended to combine the input files derived from a different Maui projects beforehand the analysis. In that way you have only to work with a single file `CalculationFileData.csv`<sup>1</sup> containing all data of the your experiment.

The herein described process provides a quick way how to combine the exported files from different Maui projects. The script `combine-sets.R`<sup>2</sup> automatically saves all combined files into the correct `input` folder. Just update the folder and subfolder names. All the rest has been taken care of for you.

1. Create in the MTXQC-project folder (e.g., `psirm_glucose/`) a new folder called `raw-data`
2. Create a subfolder for each Maui-project in `psirm_glucose/raw_data/...`
3. Copy into this folder all your Maui-derived input files altogether
4. Update the parameter of `combine-sets.R`, meaning folder name definitions, file
5. Execute the R script
6. Merged files have been generated and copied into the corresponding folder: `psirm_glucose/input-folder/gc/...` or `psirm_glucose/input-folder/inc/...`
7. Copy the renamed `ManualQuantTable.tsv` files of each Maui project into `psirm_glucose/input/quant/...`

### 3.4 Annotation-file

The annotation file relate file names with experimental conditions or specify quantification standards in your batch. Two columns - **File and Type** - are obligatory and have to be present in the annotation file. In the case of their absence `MTXQCvX_part1.Rmd` stops processing and shows an error message.

A quick way to generate an annotation file is described below:

1. Copy the first row / header of `quantMassAreaMatrix.csv` file
2. Paste & transpose the content into a new Excel-File into column A
3. Change the first entry: Metabolite -> File
4. Remove the entry `QuantMasses` at the very end of the column A
5. Add the column `Type` and specify each file either as **sample** or **addQ1\_dilution**<sup>3</sup>
6. Add more columns specifying your experimental conditions, e.g., `Cellline` and `Treatment`<sup>4</sup>
7. Save the content as `csv`-file in the `psirm_glucose/input/...`

### 3.5 Sample\_extracts-file

The `sample_extracts.csv` file is required in order to determine automatically absolute quantities in the manner of `pmol/1e+6 cells` or `pmol/mg tissue` in the `CalculationFileData.csv`.

This file requires two obligatory columns - **Extract\_vol** and **Unit**<sup>5</sup>. Please specify for each experimental condition the amount of extracted cells (count), tissue (mg) or volume of blood/plasme (ul) in the unit shown in the brackets.

The names of the columns of the experimental conditions have to match up with the annotation file. Save the file in the folder `psirm_glucose/input/....`

If the defined experimental conditions do not match up with the annotation `MTXQCvX2_part1.Rmd` exit data processing. A template file is saved for review and usage at `inst/template_files/...`

<sup>1</sup>stored in `psirm_glucose/output/quant/...`

<sup>2</sup>`inst/template_files/...`

<sup>3</sup>see for further details `additionalQuant`

<sup>4</sup>optimal: two-three parameter, max: four parameter. Consider possible combinations, e.g., `HCT116-control`, `HCT116-BPTES`

<sup>5</sup>Define: count, mg or ul

## 3.6 Internal Standard

MTXQCvX2 allows the specification of project-specific internal extraction standards. The only thing you need to do is to define the corresponding compounds as an internal standard in the `conversion_metabolite.csv` file. To do so, add `InternalStandard` to the compound in last column `Standard`.

For an classical pSIRM experiment in the Kempa lab we are using cinnamic acid. The evaluation of this compound has been integrated in Maui. Peak areas of cinnamic acid are exported from a distinct container called `cinAcid`. The exported file has to be renamed to `InternalStandard.csv` though and moved to `psirm_glucose/input/gc/....`.

If you have used a different compound as an internal extraction standard you might need to extract the peak areas of this compound from the Maui export file `quantPeakAreasMatrix.csv` file and save it in the folder `psirm_glucose/input/gc/InternalStandard.csv`, respectively. Prerequisite - you have annotated the compound in Maui.

The report of `MTXQCvX2_part1.Rmd` includes the defined internal standard for each project in a message.



## Chapter 4

# Workflow for Metmax-extracted projects

### 4.1 You want to follow this ...

- in case you have measured samples and quantification standards by GC-MS
- performed the annotation of intermediates in ChromaToF or vendor software
- exported all information into .txt files
- used metmax to extract peak areas / mass isotopomer distributions (MIDs)

### 4.2 Introduction

This document describes how to use MTXQCvX2 in combination with metmax<sup>1</sup>.

Historically, MTXQCvX2 has been developed and optimized for Maui-derived input files. The MTXQCvX2-part4.Rmd functions as a converter of metmax-derived files in order to create suitable input formats for MTXQCvX-part1.Rmd.

This module could also be used to convert tables derived from other programs as long as they are stick with the herein described table formats. Mandatory columns are referenced in the text for each kind of input file.

The general workflow of the NMTXQCvX2 project is briefly shown below in quick view. More detailed instructions are summarised in the following paragraphs.

For more detailed explanations about the individual input parameter for each module of MTXQCvX2 please proceed to read the documentation about the individual modules and their knitting parameter. The relation of knitting parameter, input and output files are described in each section.

### 4.3 Quick view

1. Generate input files: run MTXQC\_part4.Rmd<sup>2</sup>
2. Setup R-project and copy MTXQC-files
3. Knit with parameter: MTXQC\_init.Rmd
4. Copy input files into corresponding folders

---

<sup>1</sup><http://gmd.mpimp-golm.mpg.de/apps/metmax/default.htm>

<sup>2</sup>read here the instructions

5. Create annotation.csv and sample\_extracts.csv files<sup>3</sup>
6. Update metabolite names in `conversion_metabolite.csv`<sup>4</sup>
7. Define the internal standard and/or alkanes<sup>5</sup>
8. Knit with parameter: `MTXQC_ExperimentalSetup.Rmd`
9. Knit with parameter: `MTXQC_part1.Rmd`
10. Knit with parameter: `MTXQC_part2.Rmd`
11. If required - proceed with `MTXQC_part3.Rmd` for ManualValidation

## 4.4 Input files

If you need an introduction about how to use metmax - have a look at the separate documentation `Metmax_intro`.

The chapter ?? `MTXQCvX_part4` explains in detail how to use this module to generate suitable input files.

## 4.5 Annotation-file

The annotation file relate file names with experimental conditions or specify quantification standards in your batch. Two columns - **File and Type** - are obligatory and have to be present in the annotation file. In the case of their absence `MTXQCvX_part1.Rmd` stops processing and shows an error message.

A quick way to generate an annotation file is described below:

1. Copy all file names from a file of your choice
2. Paste & transpose the content into a new Excel-File into column A
3. Call column A -> File
4. Optional: Remove any non-file name entry in this column
5. Add the column Type and specify each file either as **sample**, **Q1\_dilution**, **addQ1\_dilution**<sup>6</sup>
6. Add more columns specifying your experimental conditions, e.g., Cellline and Treatment<sup>7</sup>
7. Save the content as csv-file in the `psirm_glucose/input/...`

## 4.6 Sample\_extracts-file

The `sample_extracts.csv` file is required in order to determine automatically absolute quantities in the manner of pmol/1e+6 cells or pmol/mg tissue in the `CalculationFileData.csv`.

This file requires two obligatory columns - **Extract\_vol** and **Unit**<sup>8</sup>. Please specify for each experimental condition the amount of extracted cells (count), tissue (mg) or volume of blood/plasme (ul) in the unit shown in the brackets.

The names of the columns of the experimental conditions have to match up with the annotation file. Save the file in the folder `psirm_glucose/input/....`

If the defined experimental conditions do not match up with the annotation `MTXQCvX2_part1.Rmd` exit data processing. A template file is saved for review and usage at `inst/template_files/...`

---

<sup>3</sup>Details further down this document

<sup>4</sup>Column: `Metabolite_manual`

<sup>5</sup>Also in `conversion_metabolite.csv`; see below paragraph Standards

<sup>6</sup>see for further details `additionalQuant`

<sup>7</sup>optimal: two-three parameter, max: four parameter. Consider possible combinations, e.g., HCT116-control, HCT116-BPTES

<sup>8</sup>Define: count, mg or ul

## 4.7 Update metabolite names in `conversion_metabolite.csv`

The file `conversion_metabolite.csv`, saved in `config_mtx/`, serves as a kind of translational table. It defines alternative version of metabolite library names that come in handy to plot data using shorter metabolite names. This file is also used to define settings and standard classifications. Detailed information for each column of the file are shown here: REF

### 4.7.1 Match your annotation with library names

Prior the analysis you need to match the names of your intermediates with the `conversion_metabolite.csv` file. You need to update or add the corresponding name for each intermediate in the column **Metabolite\_manual**.

General suggestion for naming conventions in ChromaToF: Metabolite\_Derivate, e.g., Lactic acid\_(2TMS). In case of the presence of main- (MP) and byproducts (BP) use: Metabolite\_Derivate\_MP/BP, e.g., Glucose\_(1MEOX)(5TMS)\_MP.

If you have annotated intermediates that are not included so far in this table please follow the instructions how to extend `conversion_metabolite.csv`.REF

### 4.7.2 Define your internal standards and alkanes

MTXQCvX2 allows the specification of project-specific internal standards. Corresponding compounds have to be marked as an internal standard in `conversion_metabolite.csv` by adding the tag **InternalStandard** in the column Standard.

If you check the box - InternalStandard in the parameter selection for `MTXQCvX2_part4.Rmd` the module searches in your input file for peak areas of the defined standard and extracts the information. It also generates the file `InternalStandard.csv` and stores it at `psirm_glucose/input/gc/...`

In the same way alkanes are defined in `conversion_metabolite.csv`. Each alkane has to be flag tagged with **Alk** in the column Standard. This gives you the opportunity to implement customized mixtures of alkanes in order to determine the retention index. `MTXQCvX_part4.Rmd` recognises the flag tag and generates `Alcane_intensities.csv` based on your input file containing peak areas and saves it in `psirm_glucose/input/gc/...`<sup>9</sup>.

The in-lab protocol considers nine alkanes from c10 to c36. Standard annotation includes an hashtag, e.g., #c10. If you use this annotation even Metmax would be able to determine the retention index.

---

<sup>9</sup>It should be alkane, I know, but Maui doesn't, unfortunately...





## Chapter 5

# MTXQCvX2 Universe

Back then in 2015 I started to think about a way how to evaluate my own datasets regarding some quality parameters. From a scratch of a few quality metrics evaluated in the frame of a very, very long Rscript until now many different generations of the MTXQC existed.

Out of this evolution cycle MTXQC developed to a suit of modules that can be used complimentary in order to evaluate and process GC-MS derived metabolomics data. `MTXQC_part1.Rmd` represents the heart of this little universe. Besides that there are many ways how to use MTXQCvX2. Proposed workflows are illustrated in chapter Workflow Maui and chapter Workflow Metmax. The flow diagram introduced here might help you to decide how to get started.

The following sections give a short overview about the main parameters and functions of each module, its generated graphical illustrations and tables.

### 5.1 MTXQC\_init.Rmd

### 5.2 MTXQC\_ExperimentalSetup

### 5.3 MTXQC\_part 1

### 5.4 MTXQC\_part 2 - Post-Processing

### 5.5 MTXQC\_part 3 - Manual Validation

### 5.6 MTXQC\_part 4 - Metmax integration



## Chapter 6

# Configuration of MTXQCvX2 - config\_mtx/...

Herein explained are the customizable tables of the MTXQCvX2 universe.

### 6.1 conversion\_metabolite.csv

Column.name

Description

Value

Metabolite\_manual

Manual defined metabolite name

#Alanine (2TMS)

Metabolite

Library name of the metabolite

Alanine\_(2TMS)\_BP\_RI:1097\_IDENT:B+C

Metabolite\_short

Short version of library name of the metabolite

Alanine\_(2TMS)

Lettercode

Lettercode version of metabolite name

Ala\_2TMS

Q1\_value

Checked if quant1:1 value available

x

Mass\_Pos

m/z-value corresponding to m\_inc

118

SE\_sel

Evaluation of the MIDs

x

Q\_sel

Evaluation for absolute quantification

x

nopsirm

Exclusively for absolute quantification

Standards

Defined as standard

InternalStandard, Alk

## 6.2 letter\_pathway\_complete.csv

Column.name

Description

Value

Letter\_\_Derivate

Derivate definition

Ala

Lettercode

Lettercode name of metabolite

Ala\_3TMS

Pathway

Ass.pathway

aa

Pathway.1

Ass. pathway - ordered for heatmap

5-aa

Met\_pathway

Ass. pathway - ordered for heatmap incl. Lettercode

5-aa\_Ala\_3TMS

Subs\_class

Substance class

aa

Met\_class  
 Substance class incl. Lettercode  
 aa\_Ala\_3TMS

### 6.3 quant1-values.csv

Column.name  
 Description  
 Value  
 Letter\_Derivate  
 Derivate name of metabolite  
 3PGA  
 Quant1\_v4  
 Quantity in (pmol)  
 43480  
 Quant1\_v3  
 Quantity in (pmol)  
 43480

### 6.4 incorporation\_calc.csv & mid\_backups.csv

Column.name  
 Description  
 Value  
 Metabolite  
 Library name of metabolite  
 Alanine\_(2TMS)\_BP\_RI:1097\_IDENT:B+C  
 Mass\_mz  
 m/z-value  
 116, 118  
 LI\_MID  
 Definition of mass level  
 m0, minc  
 Column.name  
 Description  
 Value  
 Metabolite

Library name of metabolite

Alanine\_ beta-\_(3TMS)\_MP\_RI:1435\_IDENT:A+D

Mass.m.z.

m/z value

188

BackupPeakArea

Peak area of Backup MID

4960

BackupMID

MID value for corresponding Mass.m.z.

0.8005

# Chapter 7

## pSIRM experiments

The application of stable isotopes provides a powerful tool to track the activity of metabolic pathways. the time-dependent and atom-specific routing along a metabolic pathway resolved how substrates like glucose or glutamine are used in order to maintain a certain phenotype and energetic homeostasis.

We developed an approach called pulsed stable isotope resolved metabolomics (pSIRM) enabling the quantitative evaluation of metabolite pool sizes and incorporation of stable isotopes, e.g.,  $^{13}\text{C}_6$ -glucose. A thoughtful setup of the experimental design including the applied substrates and careful experimental handling are prerequisites for a successful pSIRM experiment. Essential aspects are collected in the below paragraphs along with a number of useful tweaks.

### 7.1 Experimental design

An *in vitro* pSIRM experiment lasts in total up to three days starting from the cell seeding at day zero. Further along the way up to two media changes should be included until the application of stable isotopes and harvesting the cells maintaining the continuous availability of nutrients and avoiding the accumulation of waste products (Figure ??fig:psirm)). The media change four hours prior the harvest is set up in order to give cells time to recover from the mechanical stress of the media change. At the time point of harvest cells should be in a perfect happy state regarding metabolic environment and stress.

Choose carefully the *seeding density of your cells* in the first place. High confluency inducing contact inhibition of cell growth has a strong impact on several cellular processes including the uptake of nutrients. Try to aim for petri dishes with a maximum confluency of 75-80 %. A pre-experiment including different cell densities for seeding at a number of experimental conditions helps you to get a feeling for the cell growth in general and an expected output of cells at the time point of the harvest. Later one is useful to plan sample extraction and measurement subsequently.

For adherent cell cultures only: Include for each experimental condition an additional petri dish that is solely used to determine the cell count at the time point of your harvest. This additional plate ensures a correct determination of absolute quantities and might reduce variation of pool sizes in the statistical analysis<sup>1</sup>. Think carefully about control conditions and include cell culture dishes that are not labeled. These dishes function as a control for your labeling procedure and the natural abundance of isotopes.

---

<sup>1</sup>Pelleting these cells and snap-freezing might give useful additional samples for western blotting.




Figure 7.1: Experimental design of a pSIRM experiment distinguishing short and long labeling with stable isotopes.

## 7.2 Experimental procedures

There is a slight difference in the protocols for adherent and suspension cells. Please read instructions and footnotes carefully.

Short<sup>2</sup> and long term labeling procured differ only in the applied solvents during the labeling - either label buffer (LB) only or a combination of full label medium (LM) and label buffer (LB). Latter one is applied in order to remove extracellular metabolites of the media. Both label buffer and media contain the major nutrients / stable isotopes to keep the main substrates at constant supply at all times (Table 7.2). During the application of stable isotopes longer than a few minutes cells might sense the absence of further intermediates provided during standard cell culture procedure and adjust their metabolic program accordingly.

Solvent composition of Full label medium (LM) and label buffer (BF) for a pSIRM experiment labeling with <sup>13</sup>C-glucose.

Solvent

Base

Carbon.source

Supplements

Full label medium (LM)

DMEM, without glucose, glutamine, pyruvate

<sup>13</sup>C-Glc (2.5 g/L) <sup>12</sup>C-Gln (2 mM)

small molecules (inhibitor, antibiotics)

Label buffer (LB)

HEPES (5 mM), NaCl (140 mM), pH 7.4

<sup>13</sup>C-Glc (2.5 g/L) <sup>12</sup>C-Gln (2 mM)

The quality of your data later heavily relies on the exact handling of the cells and a *consistent timing* throughout the pSIRM experiment. Especially the step removing the LB and quenching the cells should be

<sup>2</sup>I would rather recommend this up to 2 min of labeling



a matter of a tenth of seconds rather than seconds. It is of great value to perform the cell harvest with a second person.

## 7.3 Protocol pSIRM

### 7.3.1 Adherent cell cultures

The herein described protocols are detailed explanations how to perform a pSIRM cell harvest for long term label application. If you want to label for less than 2 minutes omit solely *omit steps 6-8*.

Materials:

- Cell culture dishes, max. confluency 80 %
- Labeling media (LM) supplemented with substrates (5 ml / dish)<sup>3</sup>
- Label buffer (LB) supplemented with substrates (5 ml / dish)
- Ice-cold 50 % MeOH supplemented with 2 ug/ul cinnamic acid
- 2x 5 ml pipette and tips<sup>4</sup>
- Beaker
- Ice
- 15 ml falcons (chloroform resistant)
- Cell lifter
- Biological waste bin next to your bench

Procedure:

1. Pre-warm LB and LM in the water bath
2. Take a number of petri dishes (condition-wise including all biol. replicates)
3. Discard cell culture media
4. Carefully add *long term labeling* LM OR *short term labeling*: LB
5. Incubate cells on the bench or in an incubator
6. Discard LM (beaker)
7. Add immediatly 5 ml of LB
8. Rotate dish once in order to cover complete surface
9. Meanwhile 2nd person get prepared with 5 ml ice-cold MeOH
10. Discard LB into beaker and *immediatly* 2nd person quenches with ice-cold MeOH
11. Collect cell extracts using cell lifter
12. Transfer cell extracts into 15 ml falcons
13. Store falcons on ice until further processing

Repeat this procedure (step 6-10) for all dishes of a single condition first. Once MeOH is added metabolic processes are interrupted and cell extracts can be collected with the help of cell lifter without rush and subsequently transfered to 15 ml falcon and stored on ice until further processing (see chapter Cell extraction methanolic extracts).

Determine the cell count using your additional petri dishes for each condition.

### 7.3.2 Suspension cell cultures

Materials:

- Cell culture flasks
- Labeling media (LM) supplemented with substrates

---

<sup>3</sup>Not required for short term labeling

<sup>4</sup>Highly recommended, makes labeling and harvest super quick

- 5 ml pipette and tips<sup>5</sup>
- 1 ml pipette
- Beaker
- paper tissues
- Liquid nitrogen
- 15 ml falcons
- 1.5 ml eppendorf tubes
- Biological waste bin next to your bench

Procedure:

1. Pre-warm LM in the water bath
2. Determine the cell count of your cell suspension(s)
3. Take aliquots of  $10 - 15 \times 10^6$  cells and transfer into 15 ml falcon
4. Spin down cells very gently 300 g, 2 min at room temperature
5. Discard media into beaker
6. Resuspend cells gently in three-times 1 ml<sup>6</sup>
7. Incubate and keep warm
8. Fractionate cell suspension in three eppendorf tubes (3x 1 ml each)
9. Spin down quickly in top-bench centrifuge<sup>7</sup>
10. Discard media blandtly on paper tissues
11. Snap-freeze immediatly in liquid nitrogen
12. Store cells until further processing (see chapter Cell extraction suspension cells)

The important step here to be quick is the aliquotation of the cell suspension and subsequent spin down in the table centrifuge. Suspension cells are rather small, nevertheless  $3 \times 10^6$  cells per extract are a good starting point for GC-MS measurements.

## 7.4 Hints & notes

- The only way to be reproducible and fast is to team up with a second person.
- Keep timing consistently through the experiment.
- Keep substrate concentrations constant throughout the experiment in all solutions.
- Supplement one stable isotopic labeled substrate with all remaining substrates in non-labeled form.
- Think about nutrient levels in your cell culture and your experimental conditions. Maybe you want to change things to physiological levels.
- Add additional plate to each condition in order to have material for western blotting and others.
- Check carefully the confluency of your dishes and determine seeding densities for different conditions.
- In case of small molecule inhibitors: Try to avoid to solve them in DMSO - strong impact on chromatography.

---

<sup>5</sup>Highly recommended, makes labeling and harvest super quick

<sup>6</sup>To generate three replicates

<sup>7</sup>Most of the times 30 s are already enough

## Chapter 8

# Experimental Protocols

### 8.1 Metabolomics sample extraction

#### 8.1.1 Cell culture experiments

##### 8.1.1.1 Cell extracts in MeOH

Materials:

- Chloroform
- 15 ml falcon tubes
- thermo or rotary shaker

Procedure:

- Add 1 ml chloroform to 5 ml methanolic cell extract
- Incube for 60 min at 4°C on rotary or thermo shaker
- Centrifuge at max speed for 10 min at 4°C
- Collect polar and lipid phases into fresh falcons / tubes
- Dry under vacuum

In order to generate technical backups:

- Resuspend dried extracts in 600 ul 20% MeOH
- Shake at 4°C on thermo shaker for 30 min
- Split volumes into equal parts in fresh eppendorf tubes
- Dry under vacuum

Suggested cell density:  $2 - 3 \times 10^6$  cells / extract.

##### 8.1.1.2 Cell pellets

Materials:

- MCW suppl. cinnamic acid
- 15 ml falcon tubes
- thermo or rotary shaker

Procedure:

- Resuspend cell pellets in 1 ml MCW

- Incube for 60 min at 4°C on rotary or thermo shaker
- Add 0.5 ml of H<sub>2</sub>O to induce phase separation
- Centrifuge at max speed for 10 min, cold temperatures
- Collect polar and lipid phases into fresh falcons / tubes
- Dry under vacuum

In order to generate technical backups:

- Resuspend dried extracts in 600 ul 20% MeOH
- Shake at cold temperature on thermo shaker for 30 min
- Split volumes into equal parts in fresh eppendorf tubes
- Dry under vacuum

Suggested cell density:  $2 - 3e + 6$  cells / extract.

### 8.1.2 Tissue samples

Materials:

- MCW suppl. cinnamic acid
- ddH<sub>2</sub>O
- Eppendorf tubes
- Tissue lyzer / pulverizer

Procedure:

- Snap-freeze tissue samples
- Pulverize samples
- Aliquote 50 mg of tissue powder
- Add 1.5 ml of MCW (suppl. with cinnamic acid final conc. 2 ug/ul)
- Shake for 60 min on rotary shaker at 4°C
- Add 0.5 ml ddH<sub>2</sub>O for phase separation
- Centrifuge maximum speed, 10 min at 4°C
- Collect polar and lipid phases in fresh vessels
- Dry under vacuum

### 8.1.3 Blood samples

Material:

- Minivette for capillary blood samples
- MCW suppl. cinnamic acid
- ddH<sub>2</sub>O
- Eppendorf tubes

Procedure:

- Give 20 ul blood / sera directly into 1 ml MCW to avoid lumps
- In case of lumps sonicate samples
- Shake samples at 4°C for 800 rpm for 60 min
- Add 500 ul ddH<sub>2</sub>O and vortex shortly
- Spin down at 4°C at max speed for 10 min
- Aliquote polar phase into 2-3 times 500 ul in 1.5 ml tubes
- Aliquote lipid phase 2 times in 100 ul lower in 1.5 ml eppi
- Dry in SpeedVac (35°C)

## 8.2 Sample derivatisation

### 8.2.1 Standard protocol

Materials:

- Methoxamine (MeOx)
- Pyridine (open under the hood only!)
- MSTFA
- Alkane mix (c10-c36) in Hexane
- Chromacol vials and caps (big, small)

Mixtures:

- Solvent 1: 40 mg MeOx in 1 ml Pyridine
- Solvent 2: 10 ul Alkane mix in 1 ml MSTFA

Volumens of both solvents are shown for standard (small volumes) procedures.

Procedure:

- Make sure samples are completely dry (1 h speed vac)
- Add 20 ul (10 ul) of solvent 1 / sample
- Incubate on rotary shaker, 30°C, for 60 min
- Add 80 ul (25 ul) of solvent 2 / sample
- Incubate on rotary shaker, 37°C, for 90 min
- Centrifuge to spin down insoluble materials
- Prepare aliquotes three times 28 ul or two times 15 ul (small glass vials)
- Keep on room temperature until measurement (max. 10 days)

### 8.2.2 Specific protocols

#### 8.2.2.1 Andrash

#### 8.2.2.2 Fabian

#### 8.2.2.3 Henning



## Chapter 9

# Data processing - MAUI

### 9.1 Processing In ChromaToF

Create a new folder in ChromaToF Pegasus Acquired Samples and import your files. The processing of files for Maui-assisted annotation is a two step process. Therefore two data processing methods have to be set up and applied to all files.

#### 9.1.1 Resampling

Resampling is commonly applied and results into a data transformation enabling an improved detection of low abundant peaks and a reduction of noise. (Maybe include an example?)

The processing methods requires to tick **Export of . . . .** Subsequently, you are asked to define an output folder and the following paramter:

- Reduction rate: 4
- Beginning to end of the file
- .peg-files

#### 9.1.2 Combo-export (.cdf & .csv)

Re-import the generated .peg-files into a subfolder and apply the following data processing method.

Activate the box **asddasd** and define for both file types the following parameter.

.cdf-file:

- export directory
- 

.csv-file:

- export directory
-

## 9.2 Maui notes

## 9.3 Maui exports

With initiation of a project folder via `MTXQCvX2_init.Rmd` you created an `input`-folder containing three subfolders: `gc`, `inc`, `quant`. In the following all files that should be exported and copied into these folders are described in detail.

### 9.3.1 `input/gc/...`

Four input files are exported in order to assess the quality of the GC-MS performance of the run. The menu `Diagnostics` is selectable via right click on your Maui project name. Only the cinnamic acid peak areas are exported via the function `Export Quantification` with right click on the actual container.

Exported Files:

- `Alcane_intensities.csv` - Diagnostics/Export Alcane intensities
- `InternalStandard.csv` - cinAcid container, Export quantification, rename!
- `MassSum-73.csv` - Diagnostics/QC Mass Sum Export; enter: 73 for m/z 73
- `PeakDensities-Chroma.csv` - Diagnostics/ExportPeakDensities

### 9.3.2 `input/quant/...`

Only one container has to be exported and contains the peak areas of each metabolite and measurement. Keep in mind that you should have uploaded `pTop5 mass list` for the correct determination of peak areas in case of labeling with stable isotopes.

A further note - Maui performs absolute quantification and stores values in the `samplePeakGroups-QMQ` container. These quantities are determined by polynomial regression, and not linear regression.

The file `ManualQuantTable.tsv` is automatically generated by Maui during processing the absolute quantification.

Exported Files:

- `ManualQuantTable.tsv` - location: `Maui-project/export/QM-AbsoluteQuantification/...`<sup>1</sup>
- `quantMassAreasMatrix.csv` - Quantification export of the container `samplePeakGroups`

### 9.3.3 `input/inc/...`

This exports are only required in case of an experiment including the application of stable isotopes. It's this the case you should have performed two things:

1. Used the optional upload of `pTop5 mass list`
2. Go through the pSIRM workflow in Maui

Exported Files:

- `DataMatrix.csv` - Export % Label of container `pSIRM-samplesPeakGroups`
- `pSIRM_SpectraData.csv` - pSIRM Spectra Export of container `pSIRM-samplesPeakGroups`<sup>2</sup>

---

<sup>1</sup>Don't forget to rename it - e.g., `e17123cz`

<sup>2</sup>Requires the selection of `Natural_MIDs.txt`



## Chapter 10

# Data Processing - Metmax

10.1 Resampling

10.2 1D-basic

10.3 Reference search

10.4 Export for Metmax

10.5 Data extraction with Metmax

10.5.1 Peak areas

10.5.2 MIDs



# Chapter 11

## Solutions and Materials

### 11.1 Solutions

#### 11.1.1 Label buffer

Materials:

- ddH<sub>2</sub>O (500 ml)
- 140 mM NaCl (4.1 g)
- 5 mM Hepes (0.569 g)
- pH calibration 7.4

Procedure:

- Weigh the correct amounts of Hepes and NaCl
- Resolve in a glass bottle with 450 ml of water
- Stir carefully
- Check and adjust pH
- Adjust volumes to 500 ml

#### 11.1.2 MCW

Materials:

- Methanol
- Chloroform
- ddH<sub>2</sub>O
- Cinnamic acid stock in MeOH (2 mg/ml): final conc. 2 µg/ml

Procedure:

- Mix the solvents in the ratio of volumes - Methanol:Chloroform:Water – 5:2:1
- Supplement cinnamic acid stock 1:1000
- Store at -25°C

#### 11.1.3 Alkane-Mix

Materials:

- Hexane
- Alkanes: c10, c12, c15, c17, c19, c22, c28, c32, c36
- Thermo mixer
- Glass vials and caps

Procedure:

- Prepare stock solutions in hexane:
- c10 - c17 (liquid): 25 ul/ml
- c19 - c32: 20 mg/ml
- c36: two-times 15 mg/1.5 ml
- Warm up alkane stocks in thermo mixer 40°C
- Prepare a test mixture in equal amounts, e.g., 50 ul each, but use twice the volume of c36
- Mix test mixture with MSTFA: 10 ul / 1 ml MSTFA
- Check alkane profile by GC-MS
- If required: adjust volumes and re-test or create larger volume of your mixture for aliquots
- Store aliquots in glass vials, close well and store at 4°C
- For usage: gently warm up glass vials at 30 C on thermo mixer for 10 min and vortex before adding it to the MSTFA

Adjust the volumes of the alkane stocks in order to create a curve shaped distribution of all alkanes in the chromatogram: lower intensities for c10 and c32-36, slowly increasing intensities for the alkanes in between.

## 11.2 Materials

List of chemicals including company names and catalog number.

Chemical

Catalog.Nb

Company

Storage

Chloroform

132950-1L

Sigma

RT

Cinnamic acid

C80857 5g

Aldrich

RT

Hexane for HPLC > 97% (GC)

34859-1L

Sigma

RT

Methanol

1060351000

Merck

RT

Methoxyamine hydrochloride

226904-5G

Aldrich

RT

MSTFA - Methyltrimethylsilyltrifluoroacetamide

701270201

Macherey-Nagel

4°C

Pyridine

270970-100ML

Sigma

RT



# Chapter 12

## Output tables

### 12.1 Project-files

#### 12.1.1 MTXQC\_params.csv

Parameter

Description

Value

inputformat

Format of input files

maui

ann

File name of the annotation file

annotation.csv

sample\_ext

File name of the file containing cell counts or extracted tissue weight, volumes of blood

Sample\_extracts.csv

instd

Internal extraction standard present

TRUE

addQ

Additional quantification standards measured?

no

addQ\_Int

Integration of additional quantification standards in all batches of the MTXQC project

no

substr

Applied stable isotopes

glc

quant\_vol

Dried volume of polar phase of quantification standards in uL

500

backups

Preparation of technical aliquots

2

samples

Origin of samples

cell extracts

subf

subfolder or project folder

psirm\_glucose/

data

Nature of the experiment

pSIRM

quant

Version of quantification standards

Quant1\_v4

### 12.1.2 Maui\_params.csv

Parameter

Description

Value

spath

subfolder or project folder

psirm\_glucose/

matrix

Input file containing all peak areas

quantMassAreasMatrix.csv

mz

Input file containing m/z intensities

MassSum-73.csv

mid

Input file containing MIDs for each intermediate and file



pSIRM\_SpectraData.csv

inc\_data

Input file containing calculated stable isotope incorporation in the format 100 percent = 1

DataMatrix.csv

inputformat

Input format type

maui

intstd

If MTXQC project contains an internal extraction standard

TRUE

alkanes

If MTXQC project contains evaluation of alkane profile

TRUE

peakchroma

If MTXQC project contains evaluation of derivatisation efficiency

TRUE

mqf

If MTXQC project includes absolute quantification

TRUE

inc

If MTXQC project includes labeling with stable isotopes

TRUE

### 12.1.3 Metmax\_params.csv

## 12.2 output/gc/...

### 12.2.1 HM\_GC\_values.csv & qcmetric\_xy.csv

MTXQC exports a file summarising quality factors for each of the four parameter evaluating the GC performance. A summary representing the values illustrated in the heatmap are shown in table `HM_GC_values.csv`, individual exports for each metric in table `qcmetric_xy.csv`.

Column.name

Description

Value

Batch\_Id

Batch-Id

e18274ba

qc\_metric

QC metric factor corresponding with 1 - very good and 0 - very low

0.937254457

title

Class of QC metric

alkanes

Column.name

Description

Value

Batch\_Id

Batch-Id

e18274ba

qc\_metric

QC metric factor corresponding with 1 - very good and 0 - very low

0.937254457

title

Class of QC metric

alkanes

### 12.2.2 IntStandard\_normfactors.csv & IntStandard\_stats.csv

Column.name

Description

Value

File

File name

e18274ba\_17.cdf

PeakArea

Peak area of internal extraction standard

89308492

Batch\_Id

Bacth-Id

e18274ba

IntStd\_fac

Determined normalisation factor

0.837457514

IntStd\_eval

Evaluation of normalisation factor in relation to defined range plus/minus one standard deviation within

Column.name

Description

Value

File

File name

e18274ba\_17.cdf

PeakArea

Peak area of internal extraction standard

89308492

Batch\_Id

Batch-Id

e18274ba

IntStd\_fac

Normalisation factor

0.837457514

IntStd\_eval

Evaluation regarding QC

within

n\_batch

Number of annotated intermediates

54

mean\_batch

Mean value of internal standard peak area across the batch

106642415.3

sd\_batch

Standard deviation of the peak area across the batch of the internal extraction standard

36710894.86

### 12.2.3 Min\_Annotation.csv & SumArea\_stats.csv

Column.name

Description

Value

File

File name

e18274ba\_17.cdf

Batch\_Id

Extracted Batch-Id derived from file name

e18274ba

n\_area

Number of peak areas per file

101

sum\_area

Sum of all peak areas

44614610885

n\_total

Total number of entries (including NA)

107

mean\_batch

Mean value of peak areas across the batch

34409759414

sd\_batch

Standard deviation of peak areas across the batch

12913422445

area\_fac

Normalisation factor of the file

1.296568521

n\_50

Number of annotated peak areas per file corresponding to fifty percent of maximum number of annotated peaks in the batch

53.5

Column.name

Description

Value

Batch\_Id

Batch-Id

e18274ba

n\_50

Number corresponding to fifty percent of the maximum number of annotated peaks per file

53.5

**12.2.4 mz73\_data.csv**

Column.name

Description

Value

File

File name

e18274ba\_17.cdf

Batch\_Id

Batch-ID

e18274ba

mean\_73

Mean value of the sum of m/z 73 intensities per file

16314646.1

sd\_73

Standard deviation of the mean of the sum of m/z 73 intensities per file

143890119.5

n\_peaks

Number of intensities used for statistics

600

sum\_area

Sum of all m/z intensities

6777312761

ratio\_total

Ratio of sum\_area in relation to the sum of all annotated peak areas

0.002407244

**12.3 output/quant/...{o\_quant}****12.3.1 calcheck\_linearity.csv**

Column.name

Description

Value

Metabolite

Full library name of the metabolite

Alanine\_(3TMS)\_MP\_RI:1367\_IDENT:B+C

Batch\_Id

Batch-Id

e18274ba

File

File name

e18274ba\_53.cdf

QuantMasses

Defined quantification masses

110.0 133.0 114.0 100.0 188.0 190.0

PeakArea

Sum of peak areas based on defined QuantMasses

12710956

Type

Kind of sample

sample

Metabolite\_manual

Manual defined metabolite name

#Alanine (3TMS)

Metabolite\_short

Short version of metabolite name

Alanine\_(3TMS)

Lettercode

Lettercode version of metabolite name

Ala\_3TMS

Q1\_value

Quantity in the quantification standard Q1:1 if included in pmol

x

Mass\_Pos

m/z value that should be evaluated in case of the application of isotopes

190

SE\_sel

Intermediate used for evaluation of MIDs if present

x

Q\_sel

Intermediate used for absolute quantification if Q1:1 present

x

nopsirm

If exclusively used for quantification purposes only

Standards

Origin

Quantification standard

Qstd

adj\_r\_squared

Adjusted Rsquare value of linear regression of the calibration curve

0.997002799

slope

Slope of the calibration curve determined by linear regression

0.000578604

intercept

Intercept of the calibration curve determined by linear regression

-898.3400476

calc\_curve

Calibration curve available for absolute quantification

yes\_cal

islinear

Evaluation of peak area in relation to linear range of the calibration curve

linear

### 12.3.2 CalculationFileData.csv

Column.name

Class

Description

Value

Batch\_Id

AnnExp

Batch-Id extracted from file name

e18274ba

CL

AnnExp

Experimental parameter

BE(2)-C

Cond

AnnExp

Experimental parameter

Control

File

AnnExp

File name

e18274ba\_25.cdf

Standards

AnnExp

Defined as standard (InternalStandard, Alk)

Time

AnnExp

Experimental parameter

0

Type

AnnExp

Type of measurement

sample

Extract\_\_vol

AnnExtract

Defined extractes in count, mg or uL defined in Unit

3290000

Unit

AnnExtract

Defined unit for corresponding Extract\_\_vol

count

Lettercode

AnnMet

Lettercode version of metabolite name

Glyc3P

Metabolite

AnnMet

Library name

Glycerol-3-phosphate\_(4TMS)\_MP\_RI:1756\_IDENT:B+C

Metabolite\_manual

AnnMet

Manual defined metabolite name



#Glycerol-3-phosphate (4TMS)

Metabolite\_short

AnnMet

Short version of library name

Glycerol-3-phosphate\_(4TMS)

absconc

CalCurve

Absolute quantity in pmol

24061.77157

adj\_r\_squared

CalCurve

Adjusted Rsquare value of calibration curve

0.999667641

calc\_curve

CalCurve

Relation of peak area to calibration curve

yes\_cal

intercept

CalCurve

Intercept of calibration curve

234.0083706

islinear

CalCurve

Classification of peak area

linear

Origin

CalCurve

Source of Quant1:1 value

Qstd

Q1\_value

CalCurve

Quantity defined in Q1:1 standard

x

slope

CalCurve

Slope of calibration curve

0.000821689

corr\_absconc

ExtrFactor

Corrected absolute quantities by extraction factor

16041.18105

extr\_fac

ExtrFactor

Extraction factor derived from MTXQC\_params

0.666666667

Mass\_Pos

MTXQCsettings

m/z-value (m\_inc) used for determination of stable isotope incorporation

359

nopsirm

MTXQCsettings

Exclusivly used for absolute quantification

Q\_sel

MTXQCsettings

Used for absolute quantification

x

QuantMasses

MTXQCsettings

Defined quantification masses

101.0 299.0 133.0 357.0 103.0 359.0

SE\_sel

MTXQCsettings

Evaluation of MIDs if available

x

IntStd\_eval

Norm\_InternalStd

Evaluation of internal extraction standard

within

IntStd\_fac

Norm\_InternalStd

Normalisation factor based on internal extraction standard

1.090830836

area\_fac

Norm\_SumOfArea

Normalisation factor for sum of area normalisation

0.507591279

sum\_area

Norm\_SumOfArea

Sum of all peak areas of the file

17466093782

Conc\_microM

Quantities

Quantities in (micromolar) only determined for extracted volumes

NA

Conc\_pmo

Quantities

Quantities in (pmol)

4875.73892

IntStd\_Conc

Quantities

Quantities in (pmol) after normalisation with internal standard factor

14705.47084

IntStd\_Conc\_microM

Quantities

Quantities in (micromolar) after normalisation with internal standard factor: only for volumes

NA

IntStd\_Conc\_pmo

Quantities

Quantities in (pmol/unit) after normalisation with internal standard factor

4469.747976

IntStd\_sumA\_Conc

Quantities

Quantities in (pmol) after normalisation with internal standard factor and sum of area normalisation

28971.0865

IntStd\_sumA\_Conc\_pmo

Quantities

Quantities in (pmol/per unit) after normalisation with internal standard factor and sum of area normalisation

8805.801367

PeakArea

Quantities

Determined peak area based on QuantMasses

28998527

sumA\_Conc

Quantities

Quantities in (pmol) after sum of area normalisation

31602.5545

sumA\_Conc\_microM

Quantities

Quantities in (micromolar) after sum of area normalisation

NA

sumA\_Conc\_pmo

Quantities

Quantities in (pmol/per unit) after sum of area normalisation

9605.639664

### 12.3.3 HeatMap\_Quant\_pTop5.csv

Column.name

Description

Value

Lettercode

Short name of metabolite

Cit

Batch\_Id

Batch-Id

e18274ba

Metabolite

Library name of metabolite

Citric acid\_(4TMS)\_MP\_RI:1814\_IDENT:B+D

Par

Parameter

R2\_cal

Val

Value of the parameter for corresponding metabolite

0.996053496

**12.3.4 pTop5\_Calibration\_Samples\_lincheck.csv**

Column.name

Description

Value

Lettercode

Short name of metabolite

3PGA

islinear

Evaluation of peak area in relation to calibration curve if available

NaCal

Batch\_Id

Batch-Id

e18274ba

Origin

Origin of Quant1:1-value

NA

count

Number of peak areas

51

sum\_lin

Sum of peak area per islinear level

51

prop

Fraction of count to sum\_lin

1

**12.3.5 top5\_CalibrationInfo\_unique.csv**

Column.name

Description

Value

Metabolite

Library name of metabolite

Citric acid\_(4TMS)\_MP\_RI:1814\_IDENT:B+D

Lettercode

Lettercode name of metabolite

Cit  
 Batch\_Id  
 Batch-Id  
 e18274ba  
 Origin  
 Origin of quant1:1 value  
 Qstd  
 adj\_r\_squared  
 Adjusted Rsquare of calibration curve  
 0.996053496  
 intercept  
 Intercept of calibration curve  
 564.549288  
 slope  
 Slope of calibration curve  
 0.000194064  
 Frac\_calcurve  
 Fraction of peak area in linear range of calibration curve  
 10

### 12.3.6 top5\_QMQcurveInfo.csv

Column.name  
 Description  
 Value  
 Lettercode  
 Lettercode of metabolite name  
 Cit  
 Letter\_Derivate  
 Derivate name  
 Cit  
 Quant1\_v4  
 Quant1:1 value in (pmol)  
 52050  
 Metabolite  
 Library name of metabolite  
 Citric acid\_(4TMS)\_MP\_RI:1814\_IDENT:B+D

Batch\_Id

Batch-Id

e18274ba

Dilution

Dilution factor

0.2

ChromIntensities

Corresponding peak areas

45074572

Concentration

Concentration in (pmol)

10410

Origin

Origin of quantification standard

Qstd

Metabolite\_short

Short name of metabolite

Citric acid 275\_(4TMS)

adj\_r\_squared

Adjusted Rsquare of calibration curve

0.996053496

intercept

Intercept of calibration curve

564.549288

slope

Slope of calibration curve

0.000194064

max

Max. value of calibration curve

52050

min

Min. value of calibration curve

260.25

## 12.4 output/inc/...

### 12.4.1 HeatMap\_Incorporation.csv

Column.name

Description

Value

Lettercode

Lettercode name of metabolite

3PGA

Batch\_Id

Batch-Id

e18274ba

Par

Parameter

NA\_count

Val

Value of the parameter shown in heatmap

0.740740741

### 12.4.2 SE\_calculation\_NAscore.csv

Column.name

Description

Value

Lettercode

Lettercode name of metabolite

Cit

Batch\_Id

Batch-Id

e18274ba

na\_frac\_r

Class of NA-value

0

N

Number of MIDs

46

fracr\_prop



Fraction of MIDs with x NA-values in relation to total number of evaluated MIDs

0.851851852

### 12.4.3 SE\_classification.csv

Column.name

Description

Value

Lettercode

Lettercode name of metabolite

Cit

Batch\_Id

Batch-Id

e18274ba

File

File name

e18274ba\_\_26.cdf

sum\_spa

*expl*

4297

sum\_bpa

*expl*

1861

low3a\_ratio

*expl*

2.30897367

rel\_sb

*expl*

higher

val\_score

Validation score

confident

count\_score

MID evaluation

goodQ

**12.4.4 SE\_validation.csv**

Column.name

Description

Value

Lettercode

Lettercode name of metabolite

Cit

Batch\_Id

Batch-Id

e18274ba

count\_score

Evaluation of MID

goodQ

N\_count

Number of evaluated MIDs

54

sum\_files

Number of detected MIDs

54

prop

Proportion of MIDs for count\_score level

1

## Chapter 13

# Frequently Asked Questions

13.1 What are additional quantification standards

13.2 How do I extend `conversion_metabolite.csv`



# Chapter 14

## GC-MS settings

In the following paragraphs details of GC-MS settings are described in detail. The herein described settings have been optimized for cell extracts measured in split-mode 1:5 on the instrument Pegasus 4D-C GC-ToF-MS in 1D mode equipped with an autosampler Gerstel MPS.

### 14.0.1 Autosampler settings

The table 14.0.1 summarises the defined parameter for method settings of the Gerstel MPS.

Parameter of Gerstel MPS settings

Part

Parameter

Value

System

Runtime

40 min

System

GC Cool down time

3 min

System

Cryo Timeout

15 min

MPS - Liquid Injection

Syringe

10ulALX

MPS - Liquid Injection

Sandwich

No Sandwich

MPS - Liquid Injection

Inj. Volume (uL)

0.8

MPS - Liquid Injection

Air Volume below (uL)

1

MPS - Liquid Injection

Inj. Speed (uL/s)

20

MPS - Liquid Injection

Fill Volume (uL)

4

MPS - Liquid Injection

Fill Strokes for Sample

6

MPS - Liquid Injection

Fill Speed (uL/s)

1

MPS - Liquid Injection

Eject Speed (uL/s)

20

MPS - Liquid Injection

Viscosity Delay (s)

1

MPS - Liquid Injection

Pre Inj. Delay (s)

0

MPS - Liquid Injection

Post Inj. Delay (s)

0

MPS - Liquid Injection

Inj. Penetration (mm)

30

MPS - Liquid Injection

Sample Tray Type

VT98

MPS - Liquid Injection

Vial Penetration (mm)

32

MPS - Rinse

Wash1 Preclean

2

MPS - Rinse

Wash1 Postclean

3

MPS - Rinse

Wash1 Preclean

2

MPS - Rinse

Wash1 Postclean

3

MPS - Rinse

Sample (4.0 uL)

0

MPS - Rinse

Rinse Fill speed (uL/s)

1

MPS - Rinse

Viscosity Delay (s)

1

MPS - Rinse

Eject Speed (uL/s)

20

The injector itself is a temperature-regulated system providing the advantage of a focused injection of the sample. Gradient

## 14.0.2 Gas chromatography

The most important settings of the GC method are listed below. A graphical representation of the gradient in the oven is shown in figure 14.2 and corresponding values in table 14.0.2.

Flow path:

1. Inlet: Front
2. Capillary: GC Oven 50 m, 250 u int. diameter, 0.25 u film thickness, RTX-5 phase
3. Capillary: Detector 0.21 m, 250 u int. diameter, 0.25 u film thickness, RTX-5 phase
4. Detector: TOF

Additional parameter:




Figure 14.1: Settings temperature-regulated injection for Gerstel MPS






Figure 14.2: GC gradient - graphical representation. Rate in ( $^{\circ}\text{C}/\text{min}$ ), Target temperature in ( $^{\circ}\text{C}$ ), Duration in (min).

- Carrier Gas: Helium
- Transfer line Temperature ( $^{\circ}\text{C}$ ): 250

Front Inlet:

- Mode: split / splitless
- Flow: 1.2 mL/min, entire run
- Septum Purge Flow: 0 mL/min
- Temperature: Initial - 200 $^{\circ}\text{C}$ ; duration 0 min

Oven:

- Equilibration time (s): 60
- Modulator enabled: yes
- Modulator offset ( $^{\circ}\text{C}$ ): 5

GC gradient profile in numbers.

Rate

Target\_Temp

Duration

Initial

68

2

5

120

0

7

200

0

12

320

6

### 14.0.3 Mass spectrometer settings

# Bibliography

Pietzke, M., Zasada, C., Mudrich, S., and Kempa, S. (2014). Decoding the dynamics of cellular metabolism and the action of 3-bromopyruvate and 2-deoxyglucose using pulsed stable isotope-resolved metabolomics. *Cancer & metabolism*, 2(1):9.